

TGFβ RIII Polyclonal Conjugated Antibody

Catalog No: #C46732



Package Size: #C46732-AF350 100ul #C46732-AF405 100ul #C46732-AF488 100ul
 #C46732-AF555 100ul #C46732-AF594 100ul #C46732-AF647 100ul
 #C46732-AF680 100ul #C46732-AF750 100ul #C46732-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	TGFβ RIII Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	TGFβ RIII Polyclonal Antibody detects endogenous levels of TGFβ RIII
Immunogen Description	Synthesized peptide derived from TGFβ RIII at AA range: 731-780
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	transforming growth factor, beta receptor III
Accession No.	Swiss-Prot#:Q03167NCBI Gene ID:7049
Uniprot	Q03167
GeneID	7049;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	93
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only