

## CYFIP2 Polyclonal Conjugated Antibody

Catalog No: #C46735



Package Size: #C46735-AF350 100ul #C46735-AF405 100ul #C46735-AF488 100ul  
 #C46735-AF555 100ul #C46735-AF594 100ul #C46735-AF647 100ul  
 #C46735-AF680 100ul #C46735-AF750 100ul #C46735-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CYFIP2 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	CYFIP2 Polyclonal Antibody detects endogenous levels of CYFIP2
Immunogen Description	Synthesized peptide derived from CYFIP2 at AA range: 1171-1220
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cytoplasmic FMR1-interacting protein 2 (p53-inducible protein 121)
Accession No.	Swiss-Prot#:Q96F07NCBI Gene ID:26999
Uniprot	Q96F07
GeneID	26999;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	195
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only