

Cytochrome c Oxidase 7A2 Antibody

Catalog No: #45601



Package Size: #45601-1 50ul #45601-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome c Oxidase 7A2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Cytochrome c Oxidase 7A2 Antibody detects endogenous levels of total Cytochrome c Oxidase 7A2
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant full length protein of human COX7A2
Target Name	Cytochrome c Oxidase 7A2
Other Names	Cytochrome c oxidase subunit VIIa 3; mitochondrial; COX7AP2; COX7A3; COX7AL2; COX7S
Accession No.	Swiss-Prot#: P14406NCBI Gene ID:1347
Uniprot	P14406
GeneID	1347;
Calculated MW	11 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only