

IFN14 Antibody

Catalog No: #45617

Package Size: #45617-1 50ul #45617-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | IFN14 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | IFN14 Antibody detects endogenous levels of total IFN14 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human IFN14 |
| Target Name | IFN14 |
| Other Names | IFN-alpha-14B£B¬IFN-alpha-4B£B¬ |
| Accession No. | Swiss-Prot#: P01570/P01566/P05014/P01567NCBI Gene ID:3448/3446/3441/3444 |
| Uniprot | P01570 |
| GeneID | 3448; |
| Calculated MW | 22 kDa |
| Concentration | 1.0mg/mL |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Store at -20°C |

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only