

IFIT1B Antibody

Catalog No: #45627

Package Size: #45627-1 50ul #45627-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IFIT1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	IFIT1B Antibody detects endogenous levels of total IFIT1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IFIT1B
Target Name	IFIT1B
Other Names	IFIT1L
Accession No.	Swiss-Prot#: Q5T764NCBI Gene ID:439996
Uniprot	Q5T764
GeneID	439996;
Calculated MW	55 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only