

IGFL3 Antibody

Catalog No: #45628

Package Size: #45628-1 50ul #45628-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IGFL3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IGFL3
Target Name	IGFL3
Other Names	IGFL3, IGF-like family member 3, UNQ483, RPRC483
Accession No.	Swiss-Prot#: Q6UXB1NCBI Gene ID:388555
Uniprot	Q6UXB1
GeneID	388555;
Calculated MW	14 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only