

INSL5 Antibody

Catalog No: #45634

Package Size: #45634-1 50ul #45634-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	INSL5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	INSL5 Antibody detects endogenous levels of total INSL5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human INSL5
Target Name	INSL5
Other Names	INSL5, Insulin-like peptide INSL5, Insulin-like 5, Insulin-like peptide 5, PRO182, UNQ156
Accession No.	Swiss-Prot#: Q9Y5Q6NCBI Gene ID:10022
Uniprot	Q9Y5Q6
GeneID	10022;
Calculated MW	15 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only