

## MAGC3 Antibody

Catalog No: #45657

Package Size: #45657-1 50ul #45657-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MAGC3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MAGC3 Antibody detects endogenous levels of total MAGC3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MAGC3
Target Name	MAGC3
Other Names	HCA2
Accession No.	Swiss-Prot#: Q8TD91NCBI Gene ID:139081
Uniprot	Q8TD91
GeneID	139081;
Calculated MW	72 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only