

## MYST3 Antibody

Catalog No: #45662

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MYST3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYST3 Antibody detects endogenous levels of total MYST3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYST3
Target Name	MYST3
Other Names	KAT6A, MYST-3, K(lysine) acetyltransferase 6A, Zinc finger protein 220, ZNF220, MOZ, MYST3, RUNXBP2, ZC2HC6A
Accession No.	Swiss-Prot#: Q92794NCBI Gene ID:7994
Uniprot	Q92794
GeneID	7994;
Calculated MW	225 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only