

## STRN4 Antibody

Catalog No: #45683

Package Size: #45683-1 50ul #45683-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	STRN4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	STRN4 Antibody detects endogenous levels of total STRN4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human STRN4
Target Name	STRN4
Other Names	STRN4, Striatin-4, ZIN, Zinedin
Accession No.	Swiss-Prot#: Q9NRL3NCBI Gene ID:29888
Uniprot	Q9NRL3
GeneID	29888;
Calculated MW	80 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only