

MYCNOS Antibody

Catalog No: #45686

Package Size: #45686-1 50ul #45686-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MYCNOS Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MYCNOS Antibody detects endogenous levels of total MYCNOS
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYCNOS
Target Name	MYCNOS
Other Names	CYMN, NCYM
Accession No.	Swiss-Prot#: P40205NCBI Gene ID:10408
Uniprot	P40205
GeneID	10408;
Calculated MW	12 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only