

NKAP Antibody

Catalog No: #45687

Package Size: #45687-1 50ul #45687-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NKAP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	NKAP Antibody detects endogenous levels of total NKAP
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NKAP
Target Name	NKAP
Other Names	NKAP, NF-kappaB activating protein, NFkB activating protein, NF-kappa-B-activating protein
Accession No.	Swiss-Prot#: Q8N5F7NCBI Gene ID:79576
Uniprot	Q8N5F7
GeneID	79576;
Calculated MW	47 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only