

## GOLT1B Antibody

Catalog No: #45689

Package Size: #45689-1 50ul #45689-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GOLT1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GOLT1B Antibody detects endogenous levels of total GOLT1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GOLT1B
Target Name	GOLT1B
Other Names	GOLT1B, Golgi transport 1B, GOT1A, HGOT1a, CGI-141, Golgi transport 1 homolog B, YMR292W
Accession No.	Swiss-Prot#: Q9Y3E0NCBI Gene ID:51026
Uniprot	Q9Y3E0
GeneID	51026;
Calculated MW	15 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only