

ALKBH7 Antibody

Catalog No: #45690

Package Size: #45690-1 50ul #45690-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ALKBH7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ALKBH7 Antibody detects endogenous levels of total ALKBH7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ALKBH7
Target Name	ALKBH7
Other Names	ABH7, SPATA11
Accession No.	Swiss-Prot#: Q9BT30NCBI Gene ID:84266
Uniprot	Q9BT30
GeneID	84266;
Calculated MW	24 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only