

SOX17 Antibody

Catalog No: #45703

Package Size: #45703-1 50ul #45703-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SOX17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	SOX17 Antibody detects endogenous levels of total SOX17
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SOX17
Target Name	SOX17
Other Names	SOX17, Transcription factor SOX-17, VUR3
Accession No.	Swiss-Prot#: Q9H6I2NCBI Gene ID:64321
Uniprot	Q9H6I2
GeneID	64321;
Calculated MW	47 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

WB 1:500-2000;IHC 1:50-300

Images



Western Blot analysis of extracts from Jurkat cells, using Sox-17 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human liver cancer.

Note: This product is for in vitro research use only