

MRPL2 Antibody

Catalog No: #45714

Package Size: #45714-1 50ul #45714-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MRPL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MRPL2 Antibody detects endogenous levels of total MRPL2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MRPL2
Target Name	MRPL2
Other Names	MRPL2, CGI-22, L2mt, RPML14, MRP-L2, MRP-L14
Accession No.	Swiss-Prot#: Q5T653NCBI Gene ID:51069
Uniprot	Q5T653
GeneID	51069;
Calculated MW	33 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of MRPL2 expression in mouse lung (A) whole cell lysates.

Note: This product is for in vitro research use only