

## MRPL55 Antibody

Catalog No: #45723

Package Size: #45723-1 50ul #45723-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MRPL55 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	MRPL55 Antibody detects endogenous levels of total MRPL55
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MRPL55
Target Name	MRPL55
Other Names	MRPL55, AAVG5835, MRP-L55, PRO19675, L55nt, L55mt
Accession No.	Swiss-Prot#: Q7Z7F7NCBI Gene ID:128308
Uniprot	Q7Z7F7
GeneID	128308;
Calculated MW	15 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only