

## RPS17 Antibody

Catalog No: #45724



Package Size: #45724-1 50ul #45724-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RPS17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RPS17 Antibody detects endogenous levels of total RPS17
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RPS17
Target Name	RPS17
Other Names	RPS17, 40S ribosomal protein S17, DBA4, RPS17L1, S17, Ribosomal protein S17, RPS17L2
Accession No.	Swiss-Prot#: P08708NCBI Gene ID:6218
Uniprot	P08708
GeneID	6218;
Calculated MW	16 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only