

ABHD13 Antibody

Catalog No: #45741

Package Size: #45741-1 50ul #45741-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABHD13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ABHD13 Antibody detects endogenous levels of total ABHD13
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ABHD13
Target Name	ABHD13
Other Names	ABHD13, BA153I24.2, C13orf6, BEM46L1, RP11-153I24.2
Accession No.	Swiss-Prot#: Q7L211NCBI Gene ID:84945
Uniprot	Q7L211
GeneID	84945;
Calculated MW	39 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only