

## CHRNG Antibody

Catalog No: #45743



Package Size: #45743-1 50ul #45743-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CHRNG Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CHRNG Antibody detects endogenous levels of total CHRNG
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CHRNG
Target Name	CHRNG
Other Names	CHRNG, ACHRG, Acrg
Accession No.	Swiss-Prot#: P07510NCBI Gene ID:1146
Uniprot	P07510
GeneID	1146;
Calculated MW	58 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only