

ACBD5 Antibody

Catalog No: #45745

Package Size: #45745-1 50ul #45745-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACBD5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ACBD5 Antibody detects endogenous levels of total ACBD5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ACBD5
Target Name	ACBD5
Other Names	ACBD5, KIAA1996, Endozepine-related protein
Accession No.	Swiss-Prot#: Q5T8D3NCBI Gene ID:91452
Uniprot	Q5T8D3
GeneID	91452;
Calculated MW	60 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of ACBD5 expression in Jurkat (A), A549 (B), MCF7 (C) whole cell lysates.

Note: This product is for in vitro research use only