

ADAMTS12 Antibody

Catalog No: #45752



Package Size: #45752-1 50ul #45752-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADAMTS12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	ADAMTS12 Antibody detects endogenous levels of total ADAMTS12
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ADAMTS12
Target Name	ADAMTS12
Other Names	A disintegrin and metalloproteinase with thrombospondin motifs 12 ADAM-TS 12 ADAM-TS12 ADAMTS-12
Accession No.	Swiss-Prot#: P58397NCBI Gene ID:81792
Uniprot	P58397
GeneID	81792;
Calculated MW	178kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only