

## ADAMTS20 Antibody

Catalog No: #45753

Package Size: #45753-1 50ul #45753-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ADAMTS20 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	ADAMTS20 Antibody detects endogenous levels of total ADAMTS20
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ADAMTS20
Target Name	ADAMTS20
Other Names	A disintegrin and metalloproteinase with thrombospondin motifs 20 (EC:3.4.24.-) ADAM-TS 20 ADAM-TS20 ADAMTS-20
Accession No.	Swiss-Prot#: P59510NCBI Gene ID:80070
Uniprot	P59510
GeneID	80070;
Calculated MW	215kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only