

## C9orf23 Antibody

Catalog No: #45757

Package Size: #45757-1 50ul #45757-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	C9orf23 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	C9orf23 Antibody detects endogenous levels of total C9orf23
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human C9orf23
Target Name	C9orf23
Other Names	C9orf23
Accession No.	Swiss-Prot#: Q8N5L8NCBI Gene ID:138716
Uniprot	Q8N5L8
GeneID	138716;
Calculated MW	18 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only