

## AMOT Antibody

Catalog No: #45767



Package Size: #45767-1 50ul #45767-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | AMOT Antibody  |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Purification          | Immunogen affinity purified  |
| Applications          | WB   |
| Species Reactivity    | Hu Ms  |
| Specificity           | AMOT Antibody detects endogenous levels of total AMOT  |
| Immunogen Type        | Peptide  |
| Immunogen Description | A synthesized peptide derived from human AMOT  |
| Target Name           | AMOT   |
| Other Names           | AMOT, Angiotenin, Angiotenin p130 isoform, Angiotenin p80 isoform, KIAA1071                        |
| Accession No.         | Swiss-Prot#: Q4VCS5NCBI Gene ID:154796   |
| Uniprot               | Q4VCS5   |
| GeneID                | 154796;  |
| Calculated MW         | 118kDa   |
| Concentration         | 1.0mg/mL   |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at -20°C   |

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only