

## ATXN10 Antibody

Catalog No: #45777

Package Size: #45777-1 50ul #45777-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ATXN10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ATXN10 Antibody detects endogenous levels of total ATXN10
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ATXN10
Target Name	ATXN10
Other Names	ATXN10, Ataxin 10, Brain protein E46 homolog, E46L, HUMEEP, SCA10, Spinocerebellar ataxia 10, Ataxin-10
Accession No.	Swiss-Prot#: Q9UBB4NCBI Gene ID:25814
Uniprot	Q9UBB4
GeneID	25814;
Calculated MW	53kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of Ataxin 10 expression in Hela (A), mouse testis (B), rat testis (C) whole cell lysates.

---

Note: This product is for in vitro research use only