

## ATXN7 Antibody

Catalog No: #45779

Package Size: #45779-1 50ul #45779-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ATXN7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ATXN7 Antibody detects endogenous levels of total ATXN7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ATXN7
Target Name	ATXN7
Other Names	ATXN7, ADCAII, Ataxin 7, SCA7, OPCA3, Ataxin-7
Accession No.	Swiss-Prot#: O15265NCBI Gene ID:6314
Uniprot	O15265
GeneID	6314;
Calculated MW	95 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of Ataxin 7 expression in A549 (A), NIH3T3 (B), PC12 (C) whole cell lysates.

---

Note: This product is for in vitro research use only