

## ATP5G1 Antibody

Catalog No: #45784

Package Size: #45784-1 50ul #45784-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ATP5G1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ATP5G1 Antibody detects endogenous levels of total ATP5G1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ATP5G1
Target Name	ATP5G1
Other Names	ATP5G1, ATPase protein 9, ATPase subunit 9, ATPase subunit c, ATP synthase proteolipid P1, ATPase protein 9 isoform 1, ATP5A, ATPase subunit C isoform 1, ATP5G
Accession No.	Swiss-Prot#: P05496NCBI Gene ID:516
Uniprot	P05496
GeneID	516;
Calculated MW	14 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of ATP5G1 expression in HEK293T (A), Raw264.7 (B), H9C2 (C) whole cell lysates.

---

Note: This product is for in vitro research use only