

## MT-ATP8 Antibody

Catalog No: #45785

Package Size: #45785-1 50ul #45785-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MT-ATP8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MT-ATP8 Antibody detects endogenous levels of total MT-ATP8
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MT-ATP8
Target Name	MT-ATP8
Other Names	ATP8, ATPASE8, MTATP8
Accession No.	Swiss-Prot#: P03928NCBI Gene ID:4509
Uniprot	P03928
GeneID	4509;
Calculated MW	8 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of MT-ATP8 expression in HepG2 (A), mouse brain (B), rat stomach (C) whole cell lysates.

---

Note: This product is for in vitro research use only