

ATP5O Antibody

Catalog No: #45787

Package Size: #45787-1 50ul #45787-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ATP5O Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ATP5O Antibody detects endogenous levels of total ATP5O
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ATP5O
Target Name	ATP5O
Other Names	ATP5O, ATPO, OSCP
Accession No.	Swiss-Prot#: P48047NCBI Gene ID:539
Uniprot	P48047
GeneID	539;
Calculated MW	23 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only