

LY6A Conjugated Antibody

Catalog No: #C49480



Package Size: #C49480-AF350 100ul #C49480-AF405 100ul #C49480-AF488 100ul
 #C49480-AF555 100ul #C49480-AF594 100ul #C49480-AF647 100ul
 #C49480-AF680 100ul #C49480-AF750 100ul #C49480-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	LY6A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Species Reactivity	Ms
Immunogen Description	recombinant protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Ly-6A.2/Ly-6E. antibody Ly-6A.2/Ly-6E.1 antibody Ly6a antibody LY6A_MOUSE antibody Ly6al antibody Lymphocyte antigen 6A-2/6E-1 antibody Lymphocyte antigen Ly 6A.2/Ly 6E.1 antibody SCA-1 antibody Stem cell antigen 1 antibody T cell activating protein antibody T-cell-activating protein antibody TAP antibody
Accession No.	Swiss-Prot#:P05533
Uniprot	P05533
GeneID	110454;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	14/40 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only