

BAGE5 Antibody

Catalog No: #45799

Package Size: #45799-1 50ul #45799-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BAGE5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	BAGE5 Antibody detects endogenous levels of total BAGE5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BAGE5
Target Name	BAGE5
Other Names	BAG5, BAG-5, KIAA0873, Bcl-2-associated athanogene 5, BCL2-associated athanogene 5
Accession No.	Swiss-Prot#: Q9UL15NCBI Gene ID:9529
Uniprot	Q9UL15
GeneID	9529;
Calculated MW	51 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of extracts from Jurkat, K562 cells, using BAGE5 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only