

GCNT4 Antibody

Catalog No: #45802

Package Size: #45802-1 50ul #45802-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GCNT4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	GCNT4 Antibody detects endogenous levels of total GCNT4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GCNT4
Target Name	GCNT4
Other Names	GCNT4, Core2-GlcNAc-transferase 3, Core 2-branching enzyme 3, C2GNT3
Accession No.	Swiss-Prot#: Q9P109NCBI Gene ID:51301
Uniprot	Q9P109
GeneID	51301;
Calculated MW	53kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only