

B3GALT5 Antibody

Catalog No: #45804

Package Size: #45804-1 50ul #45804-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	B3GALT5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	B3GALT5 Antibody detects endogenous levels of total B3GALT5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human B3GALT5
Target Name	B3GALT5
Other Names	B3GALT5, B3T5, Beta-3-Gx-T5, Beta3GalT5, B3Gal-T5, B3galt5, B3GalT-V, Beta-1,3-GalTase 5, Beta3Gal-T5, GLCT5, B3GalTx
Accession No.	Swiss-Prot#: Q9Y2C3NCBI Gene ID:10317
Uniprot	Q9Y2C3
GeneID	10317;
Calculated MW	36 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only