

B4GALNT2 Antibody

Catalog No: #45806



Package Size: #45806-1 50ul #45806-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B4GALNT2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	B4GALNT2 Antibody detects endogenous levels of total B4GALNT2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human B4GALNT2
Target Name	B4GALNT2
Other Names	B4GALNT2, SD, Sda, B4GALGT2, B4GALT, GALGT2
Accession No.	Swiss-Prot#: Q8NHY0NCBI Gene ID:124872
Uniprot	Q8NHY0
GeneID	124872;
Calculated MW	63 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only