

## B4GALT6 Antibody

Catalog No: #45809

Package Size: #45809-1 50ul #45809-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	B4GALT6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	B4GALT6 Antibody detects endogenous levels of total B4GALT6
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human B4GALT6
Target Name	B4GALT6
Other Names	B4GALT6, Beta4GalT-VI, B4Gal-T6, Beta-1,4-GalTase 6, Beta4Gal-T6
Accession No.	Swiss-Prot#: Q9UBX8NCBI Gene ID:9331
Uniprot	Q9UBX8
GeneID	9331;
Calculated MW	45 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only