

B4GALT7 Antibody

Catalog No: #45810

Package Size: #45810-1 50ul #45810-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	B4GALT7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	B4GALT7 Antibody detects endogenous levels of total B4GALT7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human B4GALT7
Target Name	B4GALT7
Other Names	XGALT1
Accession No.	Swiss-Prot#: Q9UBV7NCBI Gene ID:11285
Uniprot	Q9UBV7
GeneID	11285;
Calculated MW	37 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only