

BRI3 Antibody

Catalog No: #45813

Package Size: #45813-1 50ul #45813-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BRI3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	BRI3 Antibody detects endogenous levels of total BRI3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BRI3
Target Name	BRI3
Other Names	BRI3, Brain protein I3, PRGR2, I3
Accession No.	Swiss-Prot#: O95415NCBI Gene ID:25798
Uniprot	O95415
GeneID	25798;
Calculated MW	14 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-2000

Images



Note: This product is for in vitro research use only