

CAMTA2 Antibody

Catalog No: #45832

Package Size: #45832-1 50ul #45832-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CAMTA2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CAMTA2 Antibody detects endogenous levels of total CAMTA2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CAMTA2
Target Name	CAMTA2
Other Names	CAMTA2, KIAA0909
Accession No.	Swiss-Prot#: O94983NCBI Gene ID:23125
Uniprot	O94983
GeneID	23125;
Calculated MW	132kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only