

PKC gamma Conjugated Antibody

Catalog No: #C49912



Package Size: #C49912-AF350 100ul #C49912-AF405 100ul #C49912-AF488 100ul

#C49912-AF555 100ul #C49912-AF594 100ul #C49912-AF647 100ul

#C49912-AF680 100ul #C49912-AF750 100ul #C49912-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PKC gamma Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Synthetic peptide within Human PKC gamma aa 350-450.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	KPCG_HUMAN antibody MGC57564 antibody OTTHUMP00000067291 antibody PKC-gamma antibody PKCC antibody PKCG antibody PRKCG antibody Protein kinase C gamma antibody Protein kinase C gamma polypeptide antibody Protein kinase C gamma type antibody Protein kinase C, gamma antibody SCA 14 antibody SCA14 antibody
Accession No.	Swiss-Prot#:P05129
Uniprot	P05129
GeneID	5582;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	78 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only