

CHST11 Antibody

Catalog No: #45841

Package Size: #45841-1 50ul #45841-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CHST11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CHST11 Antibody detects endogenous levels of total CHST11
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CHST11
Target Name	CHST11
Other Names	CHST11, C4ST, C4ST-1, HSA269537, IgH/CHST11 fusion, C4S-1, C4ST1
Accession No.	Swiss-Prot#: Q9NPF2NCBI Gene ID:50515
Uniprot	Q9NPF2
GeneID	50515;
Calculated MW	42 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only