

## CHST14 Antibody

Catalog No: #45845

Package Size: #45845-1 50ul #45845-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CHST14 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CHST14 Antibody detects endogenous levels of total CHST14
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CHST14
Target Name	CHST14
Other Names	CHST14, ATCS, Dermatan 4 sulfotransferase 1, D4ST-1, D4ST1, Dermatan 4-sulfotransferase 1, HD4ST1, HD4ST, HNK1ST
Accession No.	Swiss-Prot#: Q8NCH0NCBI Gene ID:113189
Uniprot	Q8NCH0
GeneID	113189;
Calculated MW	43 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

---

Note: This product is for in vitro research use only