

CA7 Antibody

Catalog No: #45846

Package Size: #45846-1 50ul #45846-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CA7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CA7 Antibody detects endogenous levels of total CA7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CA7
Target Name	CA7
Other Names	CA7, Carbonic anhydrase VII, Carbonic dehydratase VII, CAVII, Carbonate dehydratase VII, Carbonic anhydrase 7, CA-VII
Accession No.	Swiss-Prot#: P43166NCBI Gene ID:766
Uniprot	P43166
GeneID	766;
Calculated MW	30kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only