

CEACAM20 Antibody

Catalog No: #45856

Package Size: #45856-1 50ul #45856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEACAM20 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CEACAM20 Antibody detects endogenous levels of total CEACAM20
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CEACAM20
Target Name	CEACAM20
Other Names	CEACAM20, GPAD9366, UNQ9366
Accession No.	Swiss-Prot#: Q6UY09NCBI Gene ID:125931
Uniprot	Q6UY09
GeneID	125931;
Calculated MW	65 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only