

SLC7A3 Antibody

Catalog No: #45861



Package Size: #45861-1 50ul #45861-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC7A3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	SLC7A3 Antibody detects endogenous levels of total SLC7A3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SLC7A3
Target Name	SLC7A3
Other Names	SLC7A3, ATRC3, CAT3, CAT-3
Accession No.	Swiss-Prot#: Q8WY07NCBI Gene ID:84889
Uniprot	Q8WY07
GeneID	84889;
Calculated MW	67 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only