

RIC8 antibody

Catalog No: #23149

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RIC8 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 478 and 537 of Human RIC8A
Target Name	RIC8
Other Names	RIC8; MGC104517; MGC131931; MGC148073; MGC148074; synembryn
Accession No.	Swiss-Prot:Q9NPQ8Gene ID:60626
Uniprot	Q9NPQ8
GeneID	60626;
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 60kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images




Sample (30 ug of whole cell lysate)

A: A431

B: Hela

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Adrenal cortical tumor, using RIC8A antibody at 1: 100 dilution.

Note: This product is for in vitro research use only