

Phosphatase inhibitor library

Catalog No: #L9100



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------|--|
| Product Name | Phosphatase inhibitor library |
| Storage | Powder or pre-dissolved DMSO solutions in 96 well plate with optional 2D barcode Shipped with dry ice. |

Application Details

Number of Compounds: 24

Product Description

The compound library can be highly customized. You can select compounds, quantities, format (dry/solid or DMSO), plate map, and concentration to meet your specific requirement. Please check our Compound Library Instructions for more information

Note: This product is for in vitro research use only