

# TNFRSF1A Polyclonal Antibody

Catalog No: #27223



Package Size: #27223-1 50ul #27223-2 100ul


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name TNFRSF1A Polyclonal Antibody

## Images



Immunofluorescence analysis of HeLa cells using TNFRSF1A . Blue: DAPI for nuclear staining.

Note: This product is for in vitro research use only