

IKK beta Polyclonal Antibody

Catalog No: #27236

Package Size: #27236-1 50ul #27236-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IKK beta Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human IKK α ½ α ½ (NP_001547.1).
Other Names	IKKB;IKK-beta;IKK2;IKKB;IMD15;NFKBIKB;IKK beta
Accession No.	Uniprot:O14920GeneID:3551
Uniprot	O14920
GeneID	3551
Calculated MW	87kDa
SDS-PAGE MW	87kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using IKK α ½ α ½ antibody.

Background

The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.

Note: This product is for in vitro research use only