

MYLK4 Polyclonal Antibody

Catalog No: #27248

Package Size: #27248-1 50ul #27248-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MYLK4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MYLK4 (NP_001012418.2).
Other Names	MYLK4; SgK085; myosin light chain kinase family member 4
Accession No.	Swiss-Prot#:Q86YV6NCBI Gene ID:340156
Uniprot	Q86YV6
GeneID	340156;
Calculated MW	40kDa/44kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MYLK4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only